

4× EZscript Reverse Transcription Mix II (with gDNA Remover)

Catalog No.: EZB-RT2G

Description

The **EZBioscience® 4× EZscript Reverse Transcription Mix II (with gDNA Remover)** is a new-generation reverse transcription kit with higher reverse transcription efficiency compared to previous generations.

The kit is suitable for real-time RT-PCR (RT-qPCR) that contains a gDNA Remover which can effectively eliminate the contamination of genomic DNA (gDNA) or other double stranded DNA during the analysis of gene expression. In order to accurately analyze gene expression, it is necessary to detect cDNA in samples without contaminating DNA. To avoid amplification of gDNA, primers can be designed on different exons spanning introns. However, there may be cases where a suitable primer cannot be designed, as with a gene with a single exon or a gene without a long intron. Also, it may be difficult to avoid unexpected amplification from gDNA due to non-specific amplification or the existence of pseudo-genes. Moreover, some labs are heavily contaminated by the PCR products of previous tests. The kit offers a potent gDNA Remover that can eliminate double stranded DNA in RNA sample in 5 minutes at room temperature without loss of RNA. Then the first strand of cDNA is synthesized by adding the 4× EZscript RT Mix II and primers. Reaction products are applicable to subsequent PCR, qPCR and PCR cloning. The template could be mRNA, microRNA, lncRNA, circRNA, etc.

The reverse transcriptase in this Mix is a genetic engineered enzyme based on M-MLV (RNase H-) reverse transcriptase. The reverse transcriptase lacking RNase H activity is suitable for preparing full-length cDNA. And the multiple site-mutations of reverse transcriptase can obviously increase its affinity to RNA templates and its strand extending ability, which make reverse transcription reaction more efficient. Moreover, this transcriptase is rather resistant to common reverse transcriptase inhibitors. At the same time, this kit uses the latest optimized reaction system to further improve the reverse transcription efficiency. This product is also very suitable for reverse transcription using plant RNA.

Components

Components	EZB-RT2G (100 Rxns)	EZB-RT2G-L (500 Rxns)
4× EZscript RT Mix II	550 µl	550 µl × 5 tubes
Oligo dT18	110 µl	110 µl × 5 tubes
Random Hexamer	110 µl	110 µl × 5 tubes
gDNA Remover	220 µl	220 µl × 5 tubes
Nuclease free ddH ₂ O	1 ml	1 ml × 5 tubes

Storage

Store at -20°C.

Caution

Avoid RNase contamination

Please keep the environment of experiment clean. Clean gloves and mask should be worn during the experiment. Centrifuge tubes, tips and other supplies used in the experiment must be RNase-free.

Protocol

gDNA Remover treatment of RNA

1. Add 1 µg of total RNA (10 pg ~ 2 µg) or 50 ng Poly(A)+ RNA (10 pg ~ 500 ng) to a new RNase-free centrifuge tube. Add 2 µl gDNA Remover to the RNA, pipette up and down for 10 times to mix thoroughly. Incubate at room temperature (19 ~ 27°C) for 5 minutes (Optional: the RNA could be heated at 85°C for 1 min and put on ice immediately, before adding the RT reagents to it).

Reverse Transcription

2a. For microRNA, lncRNA and circRNA, Add components to the gDNA Remover-treated RNA according to the following table, then mix gently with a pipette and centrifuge the mixture briefly to the bottom of the tube.

Components	20 µl Reaction
gDNA Remover treated RNA	X µl
4× EZscript RT Mix II	5 µl
Target Specific primer (2 µM)	1 µl
Nuclease free ddH ₂ O	to 20 µl

2b. For mRNA, Add components to the gDNA Remover-treated RNA according to the following table, then, mix gently with a pipette and centrifuge the mixture briefly to the bottom of the tube (Random Hexamer is not necessary for the reaction. However, to get better qPCR results, it is recommended to add the Random Hexamer).

Components	20 µl Reaction
gDNA Remover treated RNA	X µl
4× EZscript RT Mix II	5 µl
Oligo dT18	1 µl
Random Hexamer (optional)	1 µl
Nuclease free ddH ₂ O	to 20 µl

3. Perform the reverse transcription reaction at 42°C for 15 minutes and 95°C for 30 seconds.

The cDNA products can be used in qPCR reactions immediately, or stored at -80°C for long-term storage. Avoid repeated freeze-thaw cycles.