

EZBpress Serum/Plasma RNA Purification Kit

Catalog No.: EZB-RN2

Description

The **EZBioscience®** EZBpress Serum/Plasma RNA Purification Kit provides a simple, reliable, and rapid method for isolating high-quality total RNA (including mRNA, miRNA, lncRNA, circRNA, et.al.) from human and animal serum/plasma. Serum/plasma samples are first lysed and homogenized in a strong denaturant and phenol containing buffer, which immediately inactivates RNases to ensure isolation of intact RNA. After addition of Buffer A, the homogenate is separated into aqueous and organic phases by centrifugation. RNA partitions to the upper, aqueous phase while DNA partitions to the interphase and proteins to the lower, organic phase or the interphase. The upper, aqueous phase is collected and mixed with Supplemental Reagent and ethanol. The liquid is then passed through the RNA binding Spin Column containing a silica-based membrane to which the RNA binds. Impurities are effectively removed by subsequent washing. The purified total RNA is then eluted with Elution Buffer. The purified total RNA is suitable for use in a variety of downstream applications, including: RT-PCR, RT-qPCR, RNA-seq, Northern blotting, Microarray analysis, and so on.

Components

Components	EZB-RN2 (50 Preps)
Lysis Buffer	30 ml
Buffer A	12 ml
Supplemental Reagent	550 µl
Wash Buffer 1*	8 ml
Wash Buffer 2*	8 ml
Elution Buffer	25 ml
Spin Columns (with Collection Tubes)	50 Preps

*Before using for the first time, add 32 ml of 100% ethanol to the Wash Buffer 1, 32 ml of 100% ethanol to the Wash Buffer 2.

Storage

Store the Lysis Buffer and Buffer A at 2 ~ 8°C, protect from light. Store the Supplemental Reagent at -20°C. Store other components at room temperature (When using these buffers, be careful to avoid of contamination). Divide the Elution Buffer into small aliquots upon reception is suggested.

Protocol

Sample Lysis


- Place 50 ~ 200 µl serum/plasma in a 1.5 ml centrifuge tube. Add 500 µl of Lysis Buffer. Pipette up and down for 10 times to mix thoroughly.
- After mixed, Incubate the mixture at room temperature for 5 minutes to lysis thoroughly.
- Add 150 µl Buffer A to the lysates, mix by pipetting and hand-shaking. Incubate at room temperature for 3 ~ 5 minutes.
- Centrifuge at 15,000 × g, 4°C for 3 ~ 5 minutes. Then the mixture is separated into three phases. Transfer the upper-layer supernatant (about 300 µl) to a new RNase free 1.5 ml

centrifuge tube (be careful to avoid disturbing the middle or bottom layer).

RNA Binding

- Add 10 µl Supplemental Reagent to each volume of supernatant, **then add 1.6 volume of 100% ethanol** to each volume of the above solution (supernatant + Supplemental Reagent).
- Invert the centrifuge tube for several times or pipette up and down for 10 times to mix thoroughly, and transfer the mixture to the Spin Column. Centrifuge at 4,000 × g, 4°C for 1 minute. Pour off the liquid.

Experimental Procedure Overview



Column Washing

- Add 500 µl of Wash Buffer 1 to the column. Centrifuge at 12,000 × g, 4°C for 1 minute (**be careful to avoid of contacting the bottom of the column with the liquid when taking out of the**

column). Pour off the liquid.

8. Add 500 μ l of Wash Buffer 2 to the column. Centrifuge at $12,000 \times g$, 4°C for 1 minute.
9. Pour off the liquid and eliminate the residual liquid using towel paper. Place the empty column back on the collection tube and Centrifuge at $12,000 \times g$, 4°C for 1 minute.
10. Don't need to pour off the liquid, directly transfer the column to an RNase free 1.5 ml centrifuge tube, **open the lid and keep in the air for 2 minutes**.

RNA Elution

11. Preheat the Elution Buffer to 60°C , add 20 ~ 30 μ l to the membrane in the center of the Spin Column, and then place the EP tube with the Spin Column in a water bath or metal bath at 60°C for 3 minutes to dissolve the RNA.
12. Centrifuge at $12,000 \times g$, 4°C for 1 minute (**transfer the eluate back to the column, incubate for 5 minutes and centrifuge once more will get more RNA**).
13. Discard the column, determine the RNA concentration, do the following experiment with the purified RNA, or store the RNA at -80°C until needed.

Representative Results




Figure 1. RNA isolated from serum using **EZBioscience**® EZ-press Serum/Plasma RNA Purification Kit. This qPCR detections results shows that this kit can be well used for the purification of miRNA.

Trouble shooting

The quantity of RNA purified by this Kit is too small, or the qPCR results is not so good, the Ct value is too large, or some specific genes can not be amplified normally.

Suggestions:

- a. **Examine whether any of the reagents in the Kit is contaminated:** use a new Kit as positive control to confirm whether the reagents are contaminated (If the reagents are contaminated, the results between these two will be different, obviously). Therefore, **it is recommended to divide the**

reagents into 3 ~ 4 parts, each part in a 15 ml/50 ml centrifuge tube. Be careful when taking the reagents out of the bottles to avoid of contamination.

- b. The gene specific primers for qPCR should be divided to aliquots and stored at -20°C , to decrease the possibility of degradation or contamination.

- c. **Examine whether the experiment is carried out correctly.**

e.g.:

1. During the whole process of RNA purification by this Kit, **it must be manipulated at room temperature (unless specially required) but not on ice** (until the RNA is eluted after the centrifuge), to avoid of blocking the membrane in the Spin Column by the water-insoluble substances formed during ice cooling.
2. Before using this Kit for the first time, **add 32 ml 100% ethanol to each bottle of the Wash Buffer** and mix thoroughly.
3. Transfer the upper-layer supernatant to a new tube after phase separation. **Avoid transferring any of the middle or bottom layer into the pipette when removing the upper-layer.**
4. **Add 10 μ l Supplemental Reagent to the supernatant, then add 1.6 volume of ethanol** and mix thoroughly, then load onto the column.
5. After each column washing, the column should be taken out carefully to avoid of contacting the bottom of column with liquid.
6. The column must be dried in air for 2 minutes after the centrifugation of empty column.
7. When doing the step of RNA binding and RNA elution, **load the flow-through of centrifugation to the column and repeat the centrifugation once again will get more RNA.**